

Catalog Number

68024-05P-4-3M

HPLC Column, RP Phenyl, Wide Pore, 10um, 4.0mm ID x 50mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.



Additional Info

This Cogent RP Phenyl™ Wide Pore HPLC Column has been a premium value and used by HPLC laboratories for many years. It is double end capped (fully), with traditional bonding of a 3 carbon linkage to a terminal Phenyl Ring known as Propyl Phenyl, to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L11 Classification.

Used to separate Isomers (positional), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds in Reversed Phase HPLC, these Phenyl Columns are an excellent choice.

With a 10 µm, Particle and 300A Pore Size, it works well with higher flow rates and produces lower back pressures.

This column is often chosen as scouting column for preparative methods.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$576.05

Printed On 2024/05/19