

Catalog Number

68024-15P

HPLC Column, RP Phenyl, 10um, 4.6mm ID x 150mm Length, 100A.
Cogent Brand. 1 EA.



Additional Info

This Cogent RP Phenyl™ HPLC Column has been a premium value and used by HPLC laboratories for many years. It is double end capped (fully), with traditional bonding of a 3 carbon linkage to a terminal Phenyl Ring known as Propyl Phenyl, to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L11 Classification.

Used to separate Isomers (positional), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds in Reversed Phase HPLC, these Phenyl Columns are an excellent choice.

With a 10 µm Particle, it works well with high flow rates and produces lower back pressures.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$659.85

Printed On 2024/05/19