

## Catalog Number

68024-P30-PG1-3

**Prep Guard Column Cartridge, RP Phenyl Wide Pore, 10um, 300A.**  
30mm ID x 10mm Length. Cogent Brand. 1 EA.



### Additional Info

This Guard Column Cartridge is packed with 300Å, 10µm Cogent RP C18™ Stationary Phase. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available. It is designed to protect your RP Phenyl Wide Pore, Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

*Use this Guard Column Cartridge with a 30mm ID Prep Column Holder click [HERE](#) for more information about the holder.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$537.54**

Printed On 2024/05/19