

Catalog Number

68040-05P

HPLC Column, NP Amino, 10um, 4.6mm ID x 50mm Length, 100A.
Cogent Brand. 1 EA.



Additional Info

This Cogent NP Amino™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group bonded to world class silica. This column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods.

With a 10µm Particle, this Column works well with higher flow rates producing low back pressures compared with smaller particles. Used for polar molecules with a negative charge such as Sugars and compounds with Carboxylic Acid Groups in Normal Phase or *HILIC* modes these columns are packed and shipped in Normal Phase Solvents.

This Column is often used for scouting in Preparative Chromatography.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$576.05

Printed On 2024/05/19