

## Catalog Number

68040-P50-1003M

**Process HPLC Column, NP Amino, Wide Pore, 10um, 300A. 50mm ID x 100mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent NP Amino™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group bonded to world class, wide pore silica. This column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods with larger molecules.

With a 10 µm Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures. Used for polar molecules with a negative charge such as complex Sugars and larger compounds with Carboxylic Acid Groups in Normal Phase or *HILIC* modes these columns are packed and shipped in Normal Phase Solvents.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$12,482.55**

Printed On 2024/05/19