

## Catalog Number

68100-05P-2

**Method Development Kit. Cogent RP**, Includes 5 Cogent UHPLC columns, RP C8, RP C18, RP Cyano, RP Phenyl, and RP PFP. 2.1mm x 50mm. Cogent Brand. 1 PK.



## Additional Info

This kit is an economical way to purchase 5 Cogent HPLC columns that are often used for method development and UHPLC. Each column has a 2.1 mm internal diameter and is 50 mm long, packed with 1.8  $\mu$ m particles, 100A pore size.

This kit includes the columns listed below:

- **Cogent RP™ C8**, Click [HERE](#) for more information.
- **Cogent RP™ C18**, Click [HERE](#) for more information.
- **Cogent RP™ Cyano**, Click [HERE](#) for more information.
- **Cogent RP™ Phenyl**, Click [HERE](#) for more information.
- **Cogent RP™ PFP**, Click [HERE](#) for more information.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$8,413.39**

Printed On 2024/05/20