

Catalog Number

68104-10P-2

UHPLC Column, RP C4 Butyl, 1.8 μ m, 100A. 2.1mm ID x 100mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent Reversed Phase C4, is double end capped (*fully*), with traditional bonding of a 4 carbon ligand to the highest purity silica available.

The small, 1.8 μ m Particle Size will provide high efficiency even at high flow rates and the Peak Capacity is greater than a 3 μ m Particle. Higher back pressures than with larger particle columns are typical and a special UHPLC instrument may be required for proper results.

Used for very hydrophobic or neutral compounds in Reversed Phase, this Butyl Phase is excellent for short analyses of small molecules. Expect Less retention than with a C8 column. Ideal for new method development since these columns are lower cost to run due to reduced solvent use and wastage and increased instrument throughput.

Note: A specialized UHPLC System is required for this column and your samples should be filtered with a 0.2 μ m Syringe Filter before injection.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,466.33

Printed On 2024/05/20