

Catalog Number

68322-PG2

Guard Column Cartridges, RP PFP, 3 μ m, 100A. 4.0mm ID x 10mm Length. Cogent Brand. 5/PK.



Additional Info

These Guard Column Cartridges are packed with 3 μ m, Cogent RP PFP™ Pentafluorophenyl Stationary Phase. It is double end capped (*fully*), and bonded with a four carbon link to a benzene ring with 5 Fluorines on the highest purity 100A silica available.

Designed to protect your 3.0, 4.0 or 4.6mm ID, RP PFP™, Analytical Columns from particulate and contaminants that will adsorb to it. They are packed in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

These Guard Column Cartridges must be used with a Guard Column Holder, click [HERE](#) for more information about the Holder.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$447.95

Printed On 2024/05/20