

Catalog Number

68339-PG1

Guard Column Cartridges, RP Phenyl Hexyl, 3µm, 100A. 2.1mm ID x 10mm Length. Cogent Brand. 5/PK.



Additional Info

These Guard Column Cartridges are packed with 3µm, Cogent RP Phenyl Hexyl™ Stationary Phase. They are double end capped (*fully*), with traditional bonding of a 6 carbon silane with a terminal Phenyl Ring to the highest purity silica available.

Designed to protect your 2.1mm ID, RP Phenyl Hexyl, Analytical Columns from particulate and contaminants that will adsorb to it. They are packed in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

These Guard Column Cartridges must be used with a Guard Column Holder, click [HERE](#) for more information about the Holder.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$447.95

Printed On 2024/05/20