

## Catalog Number

68340-10P

**HPLC Column, NP Amino, 3um, 4.6mm ID x 100mm Length, 100A.**  
Cogent Brand. 1 EA.



### Additional Info

This Cogent NP Amino™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group bonded to world class silica. This column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods.

With a 3µm Particle, this column is more efficient than a 5µm column of the same phase and you will experience sharper peaks. However, due to the smaller particle size, there is an increased back pressure.

Used for polar molecules with a negative charge such as Sugars and compounds with Carboxylic Acid Groups in Normal Phase or *HILIC* modes these columns are packed and shipped in Normal Phase Solvents.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$848.36**

Printed On 2024/05/19