

## Catalog Number

68340-25P-4-3M

**HPLC Column, NP Amino, Wide Pore, 3 $\mu$ m, 4.0mm ID x 250mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.**



### Additional Info

This Cogent NP Amino™ Wide Pore Stationary Phase has been a value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group, to the highest purity silica available. This makes this column ideal as an Alternate in USP Methods that use the L8 Classification, for large molecules.

With a 3  $\mu$ m Particle and 300A Pore Size, this column performs well for polar molecules with a negative charge such as complex Sugars and larger compounds with Carboxylic Acid Groups using Normal Phase or *HILIC* modes.

*These columns are packed and shipped in Normal Phase Solvents.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,036.93**

Printed On 2024/05/19