

Catalog Number

68522-P30-250

Prep HPLC Column, RP PFP, 5um, 100A. 30mm ID x 250mm Length.
Cogent Brand. 1 EA.



Additional Info

This Cogent RP PFP™ or RP Pentafluorophenyl Stationary Phase has been a premium value used by HPLC laboratories for many years and routinely called "F5". It is double end capped (*fully*) and bonded with a four carbon link to a benzene ring with 5 Fluorines on the highest purity 100A silica available.

With a 5µm Particle, this column works well with higher flow rates in preparative HPLC making this column ideal for purification of halogenated or polar samples.

A specialized HPLC System is required for this column which is packed and shipped in Reversed Phase Solvents.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$12,149.54

Printed On 2024/05/20