

## Catalog Number

68524-05P-3M

**HPLC Column, RP Phenyl, Wide Pore, 5um, 4.6mm ID x 50mm**  
Length, 300A. Cogent Brand. 1 EA.



### Additional Info

This Cogent RP Phenyl™ Wide Pore HPLC Column has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon chain with a terminal Phenyl Ring to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L11 Classification for large molecules or in new method development.

Used to separate large Isomers (*positional*), Natural Products and Related Compounds and some Peptides in Reversed Phase HPLC.

*This column is often used for scouting preparative methods.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$733.16**

Printed On 2024/05/19