

## Catalog Number

69020-05P-3

**HPLC Column, Phenyl Hydride, 4um, 3.0mm ID x 50mm Length, 100A.**  
Cogent TYPE-C Brand. 1 EA.



### Additional Info

This Cogent Phenyl Hydride HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for analysis of compounds with aromatic rings that are difficult to retain or separate in Reversed Phase, including aromatic amines, small peptides and azides.

*This size column provides similar run times as a 4.6mm ID column, but uses much less mobile phase, with lower flow rates. No special equipment is needed for this column to save solvent use.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$899.43**

Printed On 2024/05/20