

## Catalog Number

69220-HG2

**Guard Column Starter Set, Phenyl Hydride 2.o**, Complete with 5 each 2.0mm ID x 10mm Long Cartridges with 2.2um, 120A and 1 each Guard Column Holder. Cogent TYPE-C Brand. 1 PK.



### Additional Info

This Cogent Phenyl Hydride HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for purification of compounds with aromatic rings that are difficult to retain or separate in Reversed Phase, including aromatic amines, small peptides and azides.

*Use this Starter Kit with 2.1mm ID analytical columns packed with Phenyl Hydride 2.o.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,478.80**

Printed On 2024/05/19