

## Catalog Number

69220-HG3

**Guard Column Cartridge, Phenyl Hydride 2.o**, Replacement Cartridge, 2.0mm ID x 10mm Length, 2.2 $\mu$ m, 120A. Hichrom style, in an individual black case. Cogent TYPE-C Brand. 1 EA.



### Additional Info

This Cogent Phenyl Hydride HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for purification of compounds with aromatic rings that are difficult to retain or separate in Reversed Phase, including aromatic amines, small peptides and azides.

*Use this guard column cartridge with a 2.1mm ID analytical column packed with Phenyl Hydride 2.o.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$275.19**

Printed On 2024/05/19