

## Catalog Number

70000-05P

**HPLC Column, Diamond Hydride**, 4um, 4.6mm ID x 50mm Length,  
100A. Cogent TYPE-C Brand. 1 EA.



### Additional Info

This Cogent Diamond Hydride HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for analysis of polar compounds that are difficult to retain in Reversed Phase, including polar and non-polar peptides, amino acids and other biological compounds such as metabolites.

*This size column provides good resolution and short run times and is Very Popular Column used in High Through-put Screening.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$899.43**

Printed On 2024/05/19