

## Catalog Number

70200-15P-3

**HPLC Column, Diamond Hydride 2.0, 2.2 $\mu$ m, 3.0mm ID x 150mm Length, 120A. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This Cogent Diamond Hydride 2.0 HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this smaller particle, silica hydride based stationary phase provides increased Efficiency, excellent Selectivity and faster Equilibration. Considered "Near UHPLC", this column can be used on any HPLC Instrument and still deliver improved efficiency.

Chosen for analysis of polar compounds that are difficult to retain in Reversed Phase, including polar and non-polar peptides, amino acids and other biological compounds such as metabolites.

*This size column provides similar run times as a 4.6mm ID column, but uses much less mobile phase, lower flow rates and with less back pressure. No special equipment is needed to use this column.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,622.51**

Printed On 2024/05/20