

## Catalog Number

75020-10P

**HPLC Column, HPS Phenyl, 5um, 4.6mm ID x 100mm Length, 120A.**  
Cogent Brand. 1 EA.



### Additional Info

This Cogent HPS Phenyl™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is fully end capped and is manufactured with a mono functional, 4 carbon ligand that has an attached *Phenyl* group and it is bonded to the highest purity, 120A silica available making it ideal as an Alternate in USP Methods that use the L11 Classification.

With a 5µm Particle, this Column works well with higher flow rates producing lower back pressures compared to columns with smaller particles. It is often used to separate Isomers (*positional*), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds in Reversed Phase HPLC.

*Compare this column's performance to Luna™ or Inertsil™ columns.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,010.70**

Printed On 2024/05/20