

Catalog Number

82000-GH

Guard Column Holder, 10mm long TYPE-C Guard. Accepts 1.0-4.0mm ID Cartridges. Hichrom Brand. 1 EA.



General Image for Display Purposes Only

Additional Info

For Cogent TYPE-C Guard Columns (Cartridges) See Respective Column / Stationary Phases Pages on our website.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$322.34

Printed On 2024/05/20