

## Catalog Number

85318-PG1

**Guard Column Cartridges, RP C18 ABX, 3 $\mu$ m, 100A. 2.1mm ID x 10mm Length. Cogent Brand. 5/PK.**



### Additional Info

These Guard Column Cartridges are packed with 3 $\mu$ m, Cogent RP C18 ABX™ Stationary Phase. They are specially end capped (*bulky*), with traditional bonding of an 18 carbon ligand to the highest purity silica available.

Designed to protect your 2.1mm ID, RP C18 ABX, Analytical Columns from particulate and contaminants that will adsorb to it. They are packed in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

*These Guard Column Cartridges must be used with a Guard Column Holder, click [HERE](#) for more information about the Holder.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$447.95**

Printed On 2024/05/20