

Catalog Number

95020-0AV-K

Vials, Headspace, Screw Top, Glass. Amber, 20ml, 18mm precision threads, 22.5x75.5mm outer dimensions, Rounded Bottom. MicroSolv Brand. 1000/CS.



Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass, amberized and considered low actinic. Features 18mm precision threads for proper sealing.

Compatible with CTC Combi Pal®, Perkin Elmer®, Agilent® and other Headspace autosamplers. Not suitable for older versions of TurboMatrix 110® that were manufactured before 2006.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,621.85

Printed On 2024/05/20