

Catalog Number

95020-21S-T

Caps, Headspace, Screw Top, with White Silicone Rubber / Red PTFE Pre-Slit (Star Shape) Septa, fitted in Steel, Magnetic, Silver Caps with 8mm center holes. For use with 18mm Precision Thread Headspace



Additional Info

Star Shape Pre-Slit provides a greater target area for needles and are easier to penetrate but may not reseal as well as single or non-slit septa. Recommended for SPME autosamplers including CTC Combi Pal®.

Silicone Rubber / UltraClean PTFE Septa is 55° Shore A (hardness) and is 1.5mm thick. Cap has ND18 Screw Threads.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$224.15

Printed On 2024/05/19