

Catalog Number

95020-26-S

Caps, Headspace, Screw Top, with White Silicone Rubber / Blue PTFE SPME Style Septa, fitted in Steel, Magnetic, Silver Caps with 8mm center holes. For use with 18mm Precision Thread Headspace Vials.



Additional Info

Special Shape Septa with thick outer ring and thinned center provides a good sealing with the vial and a target area easier for needles to penetrate. Recommended with CTC Combi Pal® and other autosamplers with SPME or Solid Phase Micro-Extraction injectors.

Cap has ND18 Screw Threads.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$84.60

Printed On 2024/05/19