


Catalog Number

95025-01-1S

Caps, Headspace, Crimp Top, with Dark Grey Butyl Rubber Septa, fitted in Aluminum, Silver Caps with 10mm center holes. For use with 20mm Crimp Top Headspace Vials. MicroSolv Brand. 100/PK.



Additional Info

Butyl Rubber Septa is 55° Shore A (hardness) and is 3.0mm thick. Temperature resistant from -40C to 120C. Recommended for non critical analyses or when PTFE is an analyte or must be avoided.

Low cost for budget minded laboratories.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$45.54

Printed On 2024/05/20