


## Catalog Number

95025-04-2S

**Caps, Headspace**, Crimp Top, with Translucent Blue Silicone Rubber / White PTFE Septa, fitted in Aluminum, Silver Pressure Relief Caps with 8mm center holes. For use with 20mm Crimp Top Headspace Vials.



### Additional Info

Silicone Rubber / Ultraclean PTFE with Durometer of 45 Shore A, (hardness) and 3.0mm thick. Completely laminated with PTFE on both sides of the Septa which is good for sensitive analyses. The lower Durometer means that this is a softer Septa and therefore better for easier penetrations by autosampler needles. Wide temperature range from -60° to 200°C.

Clear lacquered Caps have score lines that will break open at 3.0 +/- 0.5 Bar, for safety.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$148.90**

Printed On 2024/05/19