


Catalog Number

95025-04-3D

Caps, Headspace, Crimp Top, with Translucent Blue Silicone Rubber / Transparent PTFE Septa, fitted in Steel, Magnetic, Gold Caps with 5mm center holes. For use with 20mm Crimp Top Headspace Vials. MicroSolv



Additional Info

Silicone Rubber / Ultraclean PTFE with Durometer of 45 Shore A, (hardness) and 3.2mm thick. Completely laminated with PTFE on both sides of the Septa which is good for sensitive analyses. The lower Durometer means that this is a softer Septa and therefore better for easier penetrations by autosampler needles. Wide temperature range from -60° to 200°C.

Recommended for CTC™ and other autosamplers.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$193.05

Printed On 2024/05/19