


Catalog Number

95025-10-1S

Caps, Headspace, Crimp Top, with Blue Silicone Rubber / White PTFE, High Purity Septa, fitted in Aluminum, Silver Caps with 10mm center holes. For use with 20mm Crimp Top Headspace Vials. BASIK Brand.



Additional Info

High Purity Silicone Rubber / PTFE Septa with Durometer of 35 Shore A, (hardness) and 3.0mm thick with Ultra Low Bleed Rubber that is made with a proprietary process for ultra durability, purity and low extractables.

Excellent choice for highly sensitive analyses. The very low Durometer means that this is a softer Septa and therefore better with easier penetrations by autosampler needles.

Wide temperature range from -60° to 200°C.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$62.51

Printed On 2024/05/19