


## Catalog Number

95025-14-2S

**Caps, Headspace, Crimp Top, with White Silicone Rubber / Foil (125mil) Septa, fitted in Aluminum, Silver Pressure Relief Caps with 8mm center holes. For use with 20mm Crimp Top Headspace vials. BASIK Brand.**



### Additional Info

Silicone Rubber / Aluminum Foil (125mil thick) Silver Septa with Durometer of 50 Shore A, (hardness) and 3.0mm thick is good for high temperature analyses. No color has been added to this Septa. Often used with Perkin Elmer Instruments.

Wide temperature range from -60° to 220°C. Clear lacquered Caps have score lines that will break open at 3.0 +/- 0.5 Bar, for safety.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$196.51**

Printed On 2024/05/19