

## Catalog Number

95025-23-3D

**Caps, Headspace**, Crimp Top, with White Silicone Rubber / Blue PTFE High Purity Septa, fitted in Steel, Magnetic, Gold Caps with 5mm center holes. For use with 20mm Crimp Top Headspace Vials. BASIK Brand.



### Additional Info

Silicone Rubber / PTFE with Durometer of 55 Shore A, (hardness) and 4.5mm thick. Wide temperature range from -60° to 200°C.

Recommended for CTC™ and other autosamplers.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$191.85**

Printed On 2024/05/19