


## Catalog Number

95025-23-4D

**Caps, Headspace**, Crimp Top, with White Silicone Rubber / Blue PTFE Septa, fitted in Steel, Magnetic, Gold Caps with 8mm center holes. For use with 20mm Crimp Top Headspace Vials. MicroSolv Brand. 100/PK.



### Additional Info

Compatible with SPME, Solid Phase Micro-Extraction capable autosamplers including CTC Pal®.

Silicone Rubber / UltraClean PTFE Septa is 55° Shore A (hardness) and is 1.5mm thick. Use with SPME Vials. Click [HERE](#) to view SPME Vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$181.74**

Printed On 2024/05/19