

Catalog Number

95025-E-AG10C

Vials, Headspace, Crimp Top, Glass. Clear, 10ml, 20mm Serum Finish, 22.5x46mm outer dimensions, Flat Bottom with a Long Neck. Agilent Compatible. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass. Use with aluminum or steel crimp caps for proper sealing.

Compatible with Agilent® and other Headspace autosamplers.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$63.51

Printed On 2024/05/20