

Catalog Number

95025-PE-20WC

Vials, Headspace, Crimp Top, Glass. Clear with Write-On Patch and fill lines. 20ml, 20mm Serum Finish, 23x75.5mm outer dimensions, Beveled Top and Rounded Bottom with a Tapered Neck. Perkin Elmer



Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass. Use with aluminum or steel crimp caps for proper sealing. Write-On Patch and Fill Lines are epoxy paint and is resistant to chromatography solvents.

Compatible with Perkin Elmer®, Tekmar® and other Headspace autosamplers.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$103.69

Printed On 2024/05/20