

Catalog Number

9502C-0-0CV

Vials, Crimp Top, Glass. Clear, 2ml. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass.

Use with 11mm aluminum crimp caps. These wide opening vials have *serum finish tops* that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$41.93

Printed On 2024/05/20