

Catalog Number

9502C-00-S-M

Caps, Crimp Top, with an Open Hole. Aluminum, Silver Caps. Septa are sold separately and must be inserted manually. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. BASIK Brand.



Additional Info

Click [HERE](#) for Septa page for this cap. Choose Septa for 11mm Caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$88.04

Printed On 2024/05/20