

Catalog Number

9502C-06-LCC

Caps, Snap Top, with Clear Silicone Rubber / PTFE, ultra-pure, LCMS compatible LC Clean Septa, fitted in smooth, polyethylene Blue Caps. For use with 11mm, 12x32mm, Snap Top autosampler vials. MicroSolv



Additional Info

LC Clean™, LCMS compatible.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$47.21

Printed On 2024/05/20