


Catalog Number

9502C-90C-B

Caps, Push Top, with thinned penetration point, inserted in smooth, polyethylene, Blue Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

This cap is useful for non volatile solvents or when evaporation is not considered a problem.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$18.70

Printed On 2024/05/20