

## Catalog Number

9502C-E-0AV

**Vials, Snap Top, Glass. Amber, 2ml. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass, amberized. Amber vials are often used for light sensitive compounds and considered low actinic.

*Lowest cost in class for budget minded laboratories.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$29.18**

Printed On 2024/05/21