


Catalog Number

9502C-E-15B-M

Caps, Snap Top, with 0.04" thick, White Silicone Rubber / Natural PTFE Septa, fitted in smooth, polyethylene, Light Blue Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. BASIK Brand. 1000/CS.



Additional Info

Natural PTFE is very pure and does not have dyes that could extract from it.

Our most economical Silicone Rubber / PTFE Snap Cap. Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$203.41

Printed On 2024/05/20