

## Catalog Number

9502C-E-35R-M

**Caps, Snap Top**, with 0.04" thick, White Silicone Rubber / Natural PTFE Pre-Slit Septa, fitted in smooth, polyethylene, Red Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. BASIK Brand. 1000/CS.



### Additional Info

Pre-Slit Septa are normally recommended for Waters Alliance and Acquity instruments.

Natural Septa does not have any color pigment and will not introduce extractables.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$281.23**

Printed On 2024/05/20