

Catalog Number

9502C-E-3AP

Vial & Cap Kit includes 100 EA 2ml, Snap Top, Amber Autosampler Vials and 100 EA matching Snap Caps with White Silicone Rubber / Red PTFE Pre-Slit, Soft-Guard Septa, bonded in the Blue Caps. BASIK Brand



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 11mm Vials and Caps:

Supplied with Vials: Cat. No [9502C-E-0AV](#) and Caps: Cat No. [9502C-E-30B](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$64.03

Printed On 2024/05/20