

## Catalog Number

9502C-E-43T-S

**Caps, Crimp Top, with Soft-Guard, White Silicone / Red PTFE Pre-Slit Septa, fitted in aluminum, Silver Caps.** For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. Soft-Guard Brand. 100/PK.



### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$22.76**

Printed On 2024/05/20