

Catalog Number

9502C-WDAV

Vials, Snap Top, Glass. Amber, 2ml, with Write-On Patch and fill lines. Deactivated. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK



Additional Info

First hydrolytic, Class A, Type N51A Deactivated / Silanized, pharmaceutical grade, Borosilicate Glass, amberized. Features a wide mouth opening for easy filling. Write-On Patch and fill lines are epoxy paint for reduced smearing.

Amber vials are often used for light sensitive compounds and considered low actinic.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$78.03

Printed On 2024/05/21