

Catalog Number

9502S-01-119

Caps, Screw Top, with Orange Silicone Rubber / Ivory PTFE, low durometer Septa, fitted in smooth sided, Robo Style, soft polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials.



Additional Info

Orange-Red Silicone Rubber is softer than Natural Red Rubber and is better for your autosamplers. Ivory white is inert.

Compares to National Scientific / Thermo DP™ Screw Caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$40.09

Printed On 2024/05/20