


## Catalog Number

9502S-09-SP

**Caps, Screw Top, Solid**, with an ultra-clean, Silicone Rubber / PTFE liner, in "knurled", closed Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand 100/PK.



### Additional Info

Solid top turns 9mm screw thread autosampler vials into storage vials. The liner has a silicone rubber backing to the sample facing PTFE which allows for a much better seal for volatile compounds or for solvents than PTFE only liners.

These caps are matched to work with our MicroSolv, AQ and RSA brand vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$58.20**

Printed On 2024/05/30