

## Catalog Number

9502S-10-TEF

**Caps, Screw Top**, with White Silicone Rubber / Natural TEF Septa, fitted in smooth sided, Robo Style, soft polypropylene, Blue Caps. For use with 9-425 thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

White silicone rubber and white TEF fitted septa for a lower cost cap.

TEF is equally inert as PTFE for standard chromatography.

Compares to Thermo™ / National DP™ caps.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$36.30**

Printed On 2024/05/20