

Catalog Number

9502S-10C-BLK

Caps, Screw Top, with Silicone Rubber / PTFE, low durometer Septa, bonded in "knurled", polypropylene, Black Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

Screw Thread: 9-425mm, ND9.

Non-slit for better sealing before injection. Not recommended for Waters® Alliance autosamplers. Recommended for HPLC.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$50.00

Printed On 2024/05/20