


## Catalog Number

9502S-10M-B

**Caps, Screw Top**, with White Silicone Rubber / White PTFE Septa, fitted in "knurled", polypropylene, Blue Caps. Recommended for Mass Spec. For use with 9-425mm thread, 12x32mm autosampler vials. Compares to



### Additional Info

Made with ultra-pure silicone/PTFE septa and recommended for LCMS or GCMS.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$55.21**

Printed On 2024/05/20