


## Catalog Number

9502S-11-CB

**Caps, Screw Top**, with Natural Silicone Rubber / White PTFE Septa, bonded in "knurled", polypropylene, Blue Caps. For use with 9-425 thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

Natural Silicone Rubber is made without adding color and has the same properties as White Silicone Rubber and is a less expensive bonded Septa and Cap.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$44.88**

Printed On 2024/05/20