

## Catalog Number

9502S-11-MB-K

**Caps, Screw Top**, with Tan Silicone Rubber / White PTFE Septa, bonded in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQ Brand. 1000/CS.



### Additional Info

Ultra-pure Septa and is a good choice for LCMS, GCMS & LCMS/MS. This cap is also recommended for Waters® Acquity™ autosamplers.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$416.41**

Printed On 2024/05/20